

Hip Extension

compass 540



Equipment

- Medical device class I
- Stepless height adjustment of the footplate
- Handles for left and right
- Positioning and forearm support for upper body stabilization
- Infinitely variable training pad adjustment
- "Darth frost" protective cover with storage area
- Training panel with smart operating information on the device
- Reproducible scaled positioning aids

Technical data

- Dimensions (L/W/H): 90/99-160/168 cm
- Weight block: 80 kg (16 x 5 kg)
- Total weight: up to 191 kg
- User weight: up to 250 kg
- Energy consumption (with smart assist)
 - Standby: 0.0037 kWh
 - Operation: 0.018 kWh

Optional

- [smart assist](#)
- [Proxoforce](#) software
- Fine weight adjustment 4 x 1 kg
- Movement limitation


Documents

Instruction manual compass 540

Instructions for use for older or discontinued devices are available free of charge on request from our service department at service@physiomed-group.com. Please state the exact device designation including serial number and your address. We will send them to you.


Product images






Product colors

Frame



DB 703
Micaceous iron
ore
anthracite

Housing



darth frost

Standard upholstery



F6461454
Fire

F6461197
Anthracite

Optional equipment WITH SURCHARGE

*The colors shown may differ from the actual colors depending on the monitor and screen resolution.