

Shoulder Press/Vertical Rowing

compass 600



Equipment

- Medical device class I (without isometry), class Im (with isometry)
- [smart assist](#)
- Ergonomic movement execution with very low entry weight
- Exercises can be performed with both or one arm
- Flexible grip position for exercise variations
- Fine gradation 4 x 1 kg

Illustration including optional equipment

Technical data

- Dimensions (L/W/H): 134/154/184-201 cm
- Weight block: 80 kg (16 x 5 kg)
- Total weight: up to 388 kg
- User weight: up to 200 kg
- Standby energy consumption: 0.0040 kWh
- Energy consumption operation: 0.018 kWh/0.022 kWh* (*=full equipment (position sensors))

Optional


- Isometric measuring system
- Position sensors
- Pulse receiver uncoded
- [Proxoforce](#) software
- Analog interface (data output of the force signal to an external, analog receiver system, for combination with EMG)
- Movement limitation
- 100 kg weight block (20 x 5 kg)
- Lumbar support
- Magnetic cup holder with integrated shelf for walking aids

Documents

Operating instructions
compass 600

Instructions for use for older or discontinued devices are available free of charge on request from our service department at service@physiomed-group.com. Please state the exact device designation including serial number and your address. We will send them to you.

Product images




Product colors

Frame



RAL 7021
Black gray

Housing



RAL 9006
White
aluminum

Standard upholstery



F6461454
Fire

F6461197
Anthracite

Optional equipment WITH SURCHARGE

*The colors shown may differ from the actual colors depending on the monitor and screen resolution.