

# DEEP OSCILLATION Personal Home

Deep oscillation



## Features

- Protection class: II, type BF
- Color screen
- Single-channel depth oscillation
- Frequency range from 5-250 Hz
- Contact element made of biocompatible titanium
- Comprehensive overview of all therapy parameters
- Battery operation with powerful lithium-ion rechargeable battery



- For treatment at home
- Pre-configured therapy programmes
- Intuitive application

## Technical data

- Dimensions (L/W/H): 3.1/10/19 cm
- Weight: 0.5 kg
- Input voltage; 7.2 VDC
- Input current: 1.3 ADC
- Power supply: 1 x 7.2V Li-Ion 24 Wh battery
- Output voltage; max. 400 Vs
- Load impedance: 10 MΩ
- Output frequency: 5 - 250 Hz
- Charger:
  - Mains connection: 100 - 240 VAC
  - Mains frequency: 50 - 60 Hz
  - Current consumption: max. 0.35 A

## Studies

[To the DEEP OSCILLATION studies](#)

## Documents

[Brochure](#)

[Patient information](#)

[Flyer Lymphology](#)

## Product images

