

# PHYSIOVAC-Expert

Vacuum therapy



## Features

- Protection class 1, type BF
- Logical color coding of the stimulation current and vacuum application accessories for quick and reliable assignment of the channels and the respective polarity
- Vacuum and pulsation setting in 10 levels
- Individually adjustable lower limit of the vacuum power in pulsation mode for optimum adhesion of the electrodes in every therapy situation
- Self-regulating vacuum pump for constant negative pressure
- Vibration- and noise-optimized
- Automatic switch-off by circuit breaker in the water separator (with acoustic-optical warning signal)

## Technical data

- Dimensions (L/W/H): 35/25/10 cm
- Weight: 6 kg
- Mains connection: 100 - 240 VAC
- Mains frequency: 50 - 60 Hz
- Current consumption: 0.1 A (at 230 V) or 0.2 A (at 115 V)
- Power consumption: 60 VA
- Vacuum: 0 - 0.6 bar
- Pulsations: 0 - 60 pulsations/min, adjustable in 10 steps

## Documents

[Brochure](#)

## Product images

